



TECHNICAL SUPPORT:
877.710.6999

WHAT TO EXPECT

- Each time the Unicorder is switched on, it must complete the 10 second start up routine.
- An alert will sound if the Unicorder is not receiving a good signal.

☾ DURING YOUR SLEEP STUDY

Step 1 Complete the "Watermark Medical ARES Questionnaire".

Step 2 Turn OFF the TV, radio, room fans or other electronic appliances that may interfere with the microphone.

Step 3 Switch the Unicorder ON, an audio alert will sound signifying that the Unicorder is on. After 10 seconds an audio alert will sound indicating the start-up period is complete. You can now change positions and go to sleep.

☀ IN THE MORNING

Step 1 When you wake up, if the green light on the front of the Unicorder is illuminated, switch the Unicorder OFF. If the light is no longer illuminated, the Unicorder has automatically stopped recording. This is normal, and you do not need to turn the Unicorder off.

Step 2 Return the Unicorder to the plastic bag.

Step 3 Return the Supply Box containing the Unicorder (resealed in the plastic bag), and Questionnaire to the facility where you received the ARES Unicorder.

If the same Alert sounds more than 2 times, please call Technical Support:
1.877.710.6999

ALERT GUIDE

**ALERTS DURING THE NIGHT
 DUE TO POOR SIGNAL QUALITY**

To stop an alarm, turn your head back and forth as if you are shaking your head No.

Then reposition the Unicorder on your forehead and tighten the Cannula Slip Tube.

The alarm will sound again if the problem is not properly resolved.

* If any alert sounds more than 2 times, call Technical Support.

START-UP ERROR MESSAGES

AUDIO ALERT	REPEATS EVERY	INDICATION
'The Unicorder battery is low, it must be re-charged.'	3 seconds	Battery not charged
'Call tech support, the Unicorder is not working. Code 1'	3 seconds	Hardware Problem, study cannot be performed